

Anthocyanin Content Assay Kit - Microplate Method

Product code: 112710

Product Introduction

Anthocyanins are natural pigments that are readily soluble in polar solvents and belong to the flavonoid class of compounds. They are widely found in the roots, stems, leaves, flowers, and fruits of plants, where they produce colors ranging from red to purple and are the main pigments responsible for plant coloration.

This kit uses the pH differential method to determine anthocyanin content. At pH 1.0, anthocyanins show a maximum absorption peak at 530 nm. At pH 4.5, they are converted to a colorless chalcone form and show no absorption peak at 530 nm. Based on this property, absorbance is measured at 530 nm and 700 nm under different pH conditions.

The pH differential method reduces the effects of solution pH variation and solvent differences, helping eliminate interference from non-anthocyanin substances.

Reference assay results for a blueberry sample: A1 OD values 0.408 and 0.388; A2 OD values 0.047 and 0.046; A3 OD values 0.093 and 0.090; A4 OD values 0.048 and 0.047. Actual readings may vary depending on the instrument and test conditions.

Package Contents

Pack Size	Code	Item	Volume / Quantity
100T	112710.1	Reagent I	20 mL
100T	112710.2	Reagent II	20 mL
100T	112710.3	Extraction Solution	100 mL
100T	112710.m	Manual	1 copy

Quality Standards and Safety Instructions

Raw Material and Packaging Name	Quality Standard	Main Toxicity
Reagent I	--	--
Reagent II	--	--
Extraction Solution	--	--

Transportation and Storage

Transport	Transport with ice packs.
Storage	Store at 2-8°C.
Shelf Life	180 days.

Instructions for Use

1. Anthocyanin Extraction

1. Use an extraction solution volume to sample mass ratio of 1:5 to 1:10 (mL:g).
2. It is recommended to weigh about 0.1 g of sample and add 1 mL of extraction solution.
3. Fully homogenize and transfer to an EP tube.

4. Adjust the extraction solution volume to 1 mL, cap tightly, and perform ultrasonic extraction for 2 h.
5. Centrifuge at 8000g at room temperature for 10 min.
6. Collect the supernatant for testing.

2. Assay Procedure

1. Preheat the microplate reader for at least 30 min. Preheat Reagent I and Reagent II at 25°C (room temperature) for at least 10 min.
2. Mix 20 µL of supernatant with 180 µL of Reagent I (10-fold dilution). Incubate in a 40°C water bath for 20 min, then measure absorbance at 530 nm and 700 nm. Record these as A1 and A2.
3. Mix 20 µL of supernatant with 180 µL of Reagent II (10-fold dilution). Incubate in a 40°C water bath for 20 min, then measure absorbance at 530 nm and 700 nm. Record these as A3 and A4.
4. Calculate $\Delta A = (A1 - A2) - (A3 - A4)$.

If A1 is greater than 1, increase the dilution factor while keeping the total volume at 200 µL. For example, use 10 µL supernatant and 190 µL Reagent I for a 20-fold dilution. If A1 is less than 0.1, reduce the dilution factor while keeping the total volume unchanged. For example, use 100 µL supernatant and 100 µL Reagent I for a 2-fold dilution. Keep A1 within the 0.1-1 range to improve detection sensitivity. Adjust the volume ratio of supernatant and Reagent II in the same way. Use the actual dilution factor in the calculation formula.

3. Content Calculation

$$\text{Anthocyanin content } (\mu\text{g/g fresh weight}) = [\Delta A \times V \div (\epsilon \times d) \times M \times F \times 10^6] \div W$$

$$\text{Anthocyanin content} = 33.4 \times \Delta A \times F \div W$$

V	Extraction solution volume, 1×10^{-3} L
ϵ	Molar extinction coefficient of anthocyanin, 2.69×10^4 L/mol/cm
d	Cuvette path length, 0.5 cm
M	Relative molecular mass of anthocyanin, 449.2 g/mol
F	Dilution factor
10^6	$1 \text{ g} = 10^6 \mu\text{g}$
W	Sample dry weight, g

Precautions

- Before formal testing, select 2-3 samples with large expected differences for a preliminary test.
- Required instruments and supplies not provided: microplate reader, water bath, adjustable pipette, 96-well plate, ultrasonic cleaner, mortar, and distilled water.
- The linear range of ΔA is 0.005-1.
- This product is for scientific research use by professionals only. It must not be used for clinical diagnosis or treatment, and must not be used for food or drugs or stored in ordinary residences.
- For safety and health, wear a lab coat and disposable gloves during operation.